




DAP-1 Polyclonal Antibody

Catalog No	YP-Ab-00369
Isotype	IgG
Reactivity	Human;Mouse
Applications	IHC;IF;ELISA
Gene Name	DAP
Immunogen	The antiserum was produced against synthesized peptide derived from human DAP. AA range:11-60
Specificity	DAP-1 Polyclonal Antibody detects endogenous levels of DAP-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DAP 1;DAP 1 ALPHA;DAP 1 BETA;DAP1;DLGAP1;GKAP;hGKAP;PSD 95/SAP90 binding protein 1;DAP-1;DAP;Death-associated protein 1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Epithelium,Uterus,
Function	function:Involved in mediating interferon-gamma-induced cell death.,
Background	This gene encodes a basic, proline-rich, 15-kD protein. The protein acts as a positive mediator of programmed cell death that is induced by interferon-gamma. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, May 2014],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using DAP Antibody. The picture on the right is blocked with the synthesized peptide.