




LDOC1L Polyclonal Antibody


Catalog No	YP-Ab-00432
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;IHC;IF;ELISA
Gene Name	LDOC1L
Immunogen	The antiserum was produced against synthesized peptide derived from human LDOC1L. AA range:141-190
Specificity	LDOC1L Polyclonal Antibody detects endogenous levels of LDOC1L protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/20000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LDOC1L;Protein LDOC1L;MAR6;MART6;Leucine zipper protein down-regulated in cancer cells-like;Mammalian retrotransposon-derived protein 6
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	
Tissue Specificity	Amygdala,Brain,Trachea,
Function	similarity:Belongs to the LDOC1 family.,
Background	similarity:Belongs to the LDOC1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images



Western Blot analysis of various cells using LDOC1L Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100 (4° overnight). High-pressure and temperature Tris-EDTA, pH8.0 was used for antigen retrieval. Negative control (right) obtained from antibody was pre-absorbed by immunogen peptide.



Western blot analysis of lysates from COLO205 cells, using LDOC1L Antibody. The lane on the right is blocked with the synthesized peptide.