




PEA3 Polyclonal Antibody

Catalog No	YP-Ab-00491
Isotype	IgG
Reactivity	Human;Mouse
Applications	IHC;IF;ELISA
Gene Name	ETV4
Immunogen	The antiserum was produced against synthesized peptide derived from human ETV4. AA range:281-330
Specificity	PEA3 Polyclonal Antibody detects endogenous levels of PEA3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/20000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Adenovirus E1A enhancer-binding protein;E1AF;E1A-F;ETS translocation variant 4;PEA3;ETV4;Polyomavirus enhancer activator 3 homolog;Protein PEA3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in keratinocytes.
Function	function:Transcriptional activator that binds to the enhancer of the adenovirus E1A gene; the core-binding sequence is 5'[AC]GGA[AT]GT-3'.,similarity:Belongs to the ETS family.,similarity:Contains 1 ETS DNA-binding domain.,
Background	function:Transcriptional activator that binds to the enhancer of the adenovirus E1A gene; the core-binding sequence is 5'[AC]GGA[AT]GT-3'.,similarity:Belongs to the ETS family.,similarity:Contains 1 ETS DNA-binding domain.,
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using ETV4 Antibody. The picture on the right is blocked with the synthesized peptide.