



TUG Monoclonal Antibody

Catalog No	YP-Ab-00619
Isotype	IgG
Reactivity	Mouse
Applications	WB;ELISA
Gene Name	ASPCR1
Immunogen	Purified recombinant fragment of TUG expressed in E. Coli.
Specificity	TUG Monoclonal Antibody detects endogenous levels of TUG protein.
Formulation	Purified antibody in PBS containing 0.03% sodium azide.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TUG;ASPL;ASPS;RCC17;UBXD9;UBXN9;ASPCR1;TUG/ASPCR1;ASPCR1; Tether containing UBX domain for GLUT4;Alveolar soft part sarcoma chromosomal region candidate gene 1 protein;Alveolar soft part sarcoma locus; Renal papillary cell carcinoma protein 17;UBX domain
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis using TUG Monoclonal Antibody against NIH/3T3 cell lysate.

