




GOLT1A Polyclonal Antibody

Catalog No	YP-Ab-00700
Isotype	IgG
Reactivity	Human;Mouse
Applications	IHC;IF;ELISA
Gene Name	GOLT1A
Immunogen	The antiserum was produced against synthesized peptide derived from human GOLT1A. AA range:40-89
Specificity	GOLT1A Polyclonal Antibody detects endogenous levels of GOLT1A protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/5000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GOLT1A;Vesicle transport protein GOT1A;GOT1B;Golgi transport 1 homolog A; hGOT1b
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Multi-pass membrane protein .
Tissue Specificity	Cerebellum,
Function	function:May be involved in fusion of ER-derived transport vesicles with the Golgi complex.,similarity:Belongs to the GOT1 family.,
Background	function:May be involved in fusion of ER-derived transport vesicles with the Golgi complex.,similarity:Belongs to the GOT1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using GOLT1A Antibody. The picture on the right is blocked with the synthesized peptide.