



SFT2D3 Polyclonal Antibody

Catalog No	YP-Ab-00724
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	SFT2D3
Immunogen	The antiserum was produced against synthesized peptide derived from human SFT2D3. AA range:100-150
Specificity	SFT2D3 Polyclonal Antibody detects endogenous levels of SFT2D3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SFT2 domain containing 3;SFT2D3;Vesicle transport protein SFT2C;SFT2 domain-containing protein 3
Observed Band	20-22 kDa
Calculated Molecular Weight	22 kDa
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Embryo,
Function	function:May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.,similarity:Belongs to the SFT2 family.,
Background	function:May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.,similarity:Belongs to the SFT2 family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images