




ERM Monoclonal Antibody

Catalog No	YP-Ab-00971
Isotype	IgG
Reactivity	Human
Applications	WB;ELISA
Gene Name	ETV5
Immunogen	Purified recombinant fragment of human ERM expressed in E. Coli.
Specificity	ERM Monoclonal Antibody detects endogenous levels of ERM protein.
Formulation	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Ets-related protein ERM;ETS translocation variant 5;Ets related protein ERM;ERM;ETV5
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus.
Tissue Specificity	Ubiquitous.
Function	function: Binds to DNA sequences containing the consensus nucleotide core sequence GGAA.,mass spectrometry: PubMed:11840567,similarity: Belongs to the ETS family.,similarity: Contains 1 ETS DNA-binding domain.,tissue specificity: Ubiquitous.,
Background	function: Binds to DNA sequences containing the consensus nucleotide core sequence GGAA.,mass spectrometry: PubMed:11840567,similarity: Belongs to the ETS family.,similarity: Contains 1 ETS DNA-binding domain.,tissue specificity: Ubiquitous.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis using ERM Monoclonal Antibody against Jurkat (1), NIH/3T3 (2) and MCF-7 (3) cell lysate.