




Bovine Serum Albumin Monoclonal Antibody(Mix-mA)

Catalog No	YP-Ab-01146
Isotype	IgG
Reactivity	Cow
Applications	WB
Gene Name	
Immunogen	Purified Protein
Specificity	The antibody detects endogenous Bovine Serum Albumin protein.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Bovine Serum Albumin
Observed Band	66kD
Calculated Molecular Weight	66kD
Cell Pathway	
Tissue Specificity	
Function	
Background	Albumin is the most abundant protein in the circulatory system. In mammals albumin is synthesized initially as preproalbumin by the liver. After removal of the signal peptide, the resultant proalbumin is further processed by removal of the six-residue propeptide from the new N-terminus. The albumin released into circulation possesses a half-life of 19 days.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of Bovine Serum Albumin using Bovine Serum Albumin Monoclonal Antibody.