




TLR13 Polyclonal Antibody

Catalog No	YP-Ab-01186
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	IHC;IF
Gene Name	TLR13
Immunogen	Recombinant Protein of TLR13
Specificity	The antibody detects endogenous TLR13 protein.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1:200-500. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TLR13;Toll-like receptor 13;YM3389
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	The TLRs form a family of ligand-binding proteins that engage a variety of microbial products, that on binding, activate signaling cascades leading to the synthesis of proinflammatory molecules. TLR13 is found in mouse and participates in the innate immune response to microbial agents. It acts via MYD88 and TRAF6, leading to NF-kappa-B activation, cytokine secretion and the inflammatory response.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Immunohistochemical analysis of paraffin-embedded Mouse Spleen Tissue using TLR13 Polyclonal Antibody.