



HDAC5/9 (phospho Ser259/220) Polyclonal Antibody

Catalog No	YP-Ab-01291
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;IHC;IF;ELISA
Gene Name	HDAC5/HDAC9
Immunogen	The antiserum was produced against synthesized peptide derived from human HDAC5 around the phosphorylation site of Ser259. AA range:225-274
Specificity	Phospho-HDAC5/9 (S259/220) Polyclonal Antibody detects endogenous levels of HDAC5/9 protein only when phosphorylated at S259/220.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Antigen NY CO 9;Antigen NY-CO-9;HD5;HDAC 5;HDAC5;HDAC5_HUMAN; Histone deacetylase 5;NY CO 9;HDAC5/9;HDAC5/HDAC9;Histone deacetylase 5/9;KIAA0600;HDAC9;HDAC7;HDAC7B;HDRP;KIAA0744;MITR;Histone deacetylase 9;HD9;Histone deacetylase 7B;HD7;HD7b;Histone deacetylase-related protein;MEF2-interacting transcription rep
Observed Band	121kD
Calculated Molecular Weight	121kD
Cell Pathway	Nucleus. Cytoplasm. Shuttles between the nucleus and the cytoplasm. In muscle cells, it shuttles into the cytoplasm during myocyte differentiation. The export to cytoplasm depends on the interaction with a 14-3-3 chaperone protein and is due to its phosphorylation at Ser-259 and Ser-498 by AMPK, CaMK1 and SIK1.
Tissue Specificity	Ubiquitous.
Function	catalytic activity:Hydrolysis of an N(6)-acetyl-lysine residue of a histone to yield a deacetylated histone.,domain:The nuclear export sequence mediates the shuttling between the nucleus and the cytoplasm.,function:Responsible for the deacetylation of lysine residues on the N-terminal part of the core histones (H2A,



H2B, H3 and H4). Histone deacetylation gives a tag for epigenetic repression and plays an important role in transcriptional regulation, cell cycle progression and developmental events. Histone deacetylases act via the formation of large multiprotein complexes. Involved in muscle maturation by repressing transcription of myocyte enhancer MEF2C. During muscle differentiation, it shuttles into the cytoplasm, allowing the expression of myocyte enhancer factors.,PTM:Phosphorylated by CaMK at Ser-259 and Ser-498. The phosphorylation is required for the export to the cytoplasm.,PTM:

Background

Histones play a critical role in transcriptional regulation, cell cycle progression, and developmental events. Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by this gene belongs to the class II histone deacetylase/acuc/apha family. It possesses histone deacetylase activity and represses transcription when tethered to a promoter. It coimmunoprecipitates only with HDAC3 family member and might form multicomplex proteins. It also interacts with myocyte enhancer factor-2 (MEF2) proteins, resulting in repression of MEF2-dependent genes. This gene is thought to be associated with colon cancer. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],

matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions


This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images




Western Blot analysis of HepG2 cells using
Phospho-HDAC5/9 (S259/220) Polyclonal Antibody
diluted at 1:2000




Western Blot analysis of COLO HELA 293T cells using
Phospho-HDAC5/9 (S259/220) Polyclonal Antibody
diluted at 1:2000



Enzyme-Linked Immunosorbent Assay
(Phospho-ELISA) for Immunogen Phosphopeptide
(Phospho-left) and Non-Phosphopeptide
(Phospho-right), using HDAC5 (Phospho-Ser259)
Antibody



Western blot analysis of lysates from HepG2 cells,
using HDAC5 (Phospho-Ser259) Antibody. The lane on
the right is blocked with the phospho peptide.



Immunohistochemical analysis of paraffin-embedded
human Colon cancer. 1, Antibody was diluted at
1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for
antigen retrieval. 3, Secondary antibody was diluted at
1:200(room temperature, 45min).