



NFAT5 (phospho Ser155) Polyclonal Antibody

Catalog No	YP-Ab-01329
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;IHC
Gene Name	NFAT5
Immunogen	The antiserum was produced against synthesized peptide derived from human NFAT5/TonEBP around the phosphorylation site of Ser155. AA range:121-170
Specificity	Phospho-NFAT5 (S155) Polyclonal Antibody detects endogenous levels of NFAT5 protein only when phosphorylated at S155.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Glutamine rich protein H65;KIAA0827;NF AT5;NF-AT5;NFAT 5;NFAT L1;NFAT like protein 1;NFAT5;NFAT5_HUMAN;NFATL 1;NFATL1;NFATZ;Nuclear factor of activated T cells 5;Nuclear factor of activated T cells 5 tonicity responsive; Nuclear factor of activated T cells;Nuclear factor of activated T-cells 5;OREBP; Osmotic response element binding protein;T cell transcription factor NFAT 5;T cell transcription factor NFAT5;T-cell transcription factor NFAT5;TonE binding protein;TonE-binding protein;TonEBP;Tonicity responsive enhancer binding protein;Tonicity-responsive enhancer-binding protein
Observed Band	200kD
Calculated Molecular Weight	200kD
Cell Pathway	Nucleus . Cytoplasm . Nuclear distribution increases under hypertonic conditions. .
Tissue Specificity	Widely expressed, with highest levels in skeletal muscle, brain, heart and peripheral blood leukocytes.
Function	alternative products:Experimental confirmation may be lacking for some isoforms,function:Plays a role in the inducible expression of genes. Regulates hypertonicity-induced cellular accumulation of osmolytes.,similarity:Contains 1



RHD (Rel-like) domain.,subunit:Does not bind with Fos and Jun transcription factors. But might be capable of forming stable dimers with DNA elements., tissue specificity:Highest levels in skeletal muscle, brain, heart and peripheral blood leukocytes. Also expressed in placenta, lung, liver, kidney, pancreas, spleen, thymus, prostate, testis, ovary, small intestine and colon.,

Background

The product of this gene is a member of the nuclear factors of activated T cells family of transcription factors. Proteins belonging to this family play a central role in inducible gene transcription during the immune response. This protein regulates gene expression induced by osmotic stress in mammalian cells. Unlike monomeric members of this protein family, this protein exists as a homodimer and forms stable dimers with DNA elements. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],

matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions


This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images



Enzyme-Linked Immunosorbent Assay
(Phospho-ELISA) for Immunogen Phosphopeptide
(Phospho-left) and Non-Phosphopeptide
(Phospho-right), using NFAT5/TonEBP
(Phospho-Ser155) Antibody



Western blot analysis of lysates from HeLa cells treated with forskolin 40nM 30', using NFAT5/TonEBP (Phospho-Ser155) Antibody. The lane on the right is blocked with the phospho peptide.



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).