



HAND1 (phospho Ser98) Polyclonal Antibody

Catalog No	YP-Ab-01406
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	IHC;IF;ELISA
Gene Name	HAND1
Immunogen	The antiserum was produced against synthesized peptide derived from human HAND1 around the phosphorylation site of Ser98. AA range:71-120
Specificity	Phospho-HAND1 (S98) Polyclonal Antibody detects endogenous levels of HAND1 protein only when phosphorylated at S98.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	autonomic nervous system and neural crest derivatives-expressed protein 1; Basic helix loop helix transcription factor HAND1;bHLHa27;Class A basic helix-loop-helix protein 27;eHAND;Extraembryonic tissues;Extraembryonic tissues heart autonomic nervous system and neural crest derivatives expressed protein 1;HAND 1;HAND1;HAND1_HUMAN;Heart and neural crest derivatives expressed 1;Heart and neural crest derivatives expressed protein 1;heart;Heart-and neural crest derivatives-expressed protein 1;Hxt;Thing 1;Thing1;eH
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleoplasm . Nucleus, nucleolus . Interaction with MDFIC sequesters it into the nucleolus, preventing the transcription factor activity. Phosphorylation by PLK4 disrupts the interaction with MDFIC and releases it from the nucleolus, leading to transcription factor activity (By similarity). .
Tissue Specificity	Heart.
Function	function:Plays an essential role in early trophoblast differentiation and in cardiac morphogenesis. In the adult, could be required for ongoing expression of cardiac-specific genes. Binds the DNA sequence 5'-NRTCTG-3' (non-canonical



E-box).,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,
subunit:Efficient DNA binding requires dimerization with another bHLH protein.
Forms homodimers and heterodimers with TCF3 gene products E12 and E47,
HAND2 and HEY1, HEY2 and HEYL (hairy-related transcription factors).,tissue
specificity:Heart.,

Background

The protein encoded by this gene belongs to the basic helix-loop-helix family of transcription factors. This gene product is one of two closely related family members, the HAND proteins, which are asymmetrically expressed in the developing ventricular chambers and play an essential role in cardiac morphogenesis. Working in a complementary fashion, they function in the formation of the right ventricle and aortic arch arteries, implicating them as mediators of congenital heart disease. In addition, it has been suggested that this transcription factor may be required for early trophoblast differentiation. [provided by RefSeq, Jul 2008],



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100(4° overnight). High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtained from antibody was pre-absorbed by immunogen peptide.