




BTBD6 Polyclonal Antibody


Catalog No	YP-Ab-01576
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;IHC;IF;ELISA
Gene Name	BTBD6
Immunogen	The antiserum was produced against synthesized peptide derived from human BTBD6. AA range:91-140
Specificity	BTBD6 Polyclonal Antibody detects endogenous levels of BTBD6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BDPL;BTB (POZ) domain containing 6;BTBD6;Lens BTB domain protein;BTB/POZ domain-containing protein 6
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	Cytoplasm . Found in punctated bodies in the cytoplasm. .
Tissue Specificity	Expressed in lens.
Function	similarity:Contains 1 BTB (POZ) domain.,tissue specificity:Expressed in lens.,
Background	similarity:Contains 1 BTB (POZ) domain.,tissue specificity:Expressed in lens.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using BTBD6 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night