




# BTF3L4 Polyclonal Antibody

|                                    |                                                                                                                                   |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-Ab-01577                                                                                                                       |
| <b>Isotype</b>                     | IgG                                                                                                                               |
| <b>Reactivity</b>                  | Human;Mouse                                                                                                                       |
| <b>Applications</b>                | IHC;IF;ELISA                                                                                                                      |
| <b>Gene Name</b>                   | BTF3L4                                                                                                                            |
| <b>Immunogen</b>                   | The antiserum was produced against synthesized peptide derived from human BTF3L4. AA range:94-143                                 |
| <b>Specificity</b>                 | BTF3L4 Polyclonal Antibody detects endogenous levels of BTF3L4 protein.                                                           |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                           |
| <b>Source</b>                      | Polyclonal, Rabbit,IgG                                                                                                            |
| <b>Purification</b>                | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.             |
| <b>Dilution</b>                    | IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200                                                                                  |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                           |
| <b>Purity</b>                      | ≥90%                                                                                                                              |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                      |
| <b>Synonyms</b>                    | BTF3L4;RP4 800M22.5;Transcription factor BTF3 homolog 4;Basic transcription factor 3-like 4                                       |
| <b>Observed Band</b>               | 17kD                                                                                                                              |
| <b>Calculated Molecular Weight</b> | 17kD                                                                                                                              |
| <b>Cell Pathway</b>                |                                                                                                                                   |
| <b>Tissue Specificity</b>          | Lymph,Placenta,Prostate,                                                                                                          |
| <b>Function</b>                    | similarity:Belongs to the NAC-beta family.,similarity:Contains 1 NAC-A/B (NAC-alpha/beta) domain.,                                |
| <b>Background</b>                  | similarity:Belongs to the NAC-beta family.,similarity:Contains 1 NAC-A/B (NAC-alpha/beta) domain.,                                |
| <b>matters needing attention</b>   | Avoid repeated freezing and thawing!                                                                                              |
| <b>Usage suggestions</b>           | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



## Products Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using BTF3L4 Antibody. The picture on the right is blocked with the synthesized peptide.