



# DDX19B Polyclonal Antibody

|                                    |                                                                                                                                                                                                                                                                                                                                                              |
|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-Ab-01639                                                                                                                                                                                                                                                                                                                                                  |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                                                                          |
| <b>Reactivity</b>                  | Human;Mouse                                                                                                                                                                                                                                                                                                                                                  |
| <b>Applications</b>                | WB;IHC;IF;ELISA                                                                                                                                                                                                                                                                                                                                              |
| <b>Gene Name</b>                   | DDX19B                                                                                                                                                                                                                                                                                                                                                       |
| <b>Immunogen</b>                   | The antiserum was produced against synthesized peptide derived from human DDX19B. AA range:1-50                                                                                                                                                                                                                                                              |
| <b>Specificity</b>                 | DDX19B Polyclonal Antibody detects endogenous levels of DDX19B protein.                                                                                                                                                                                                                                                                                      |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                                                                                                      |
| <b>Source</b>                      | Polyclonal, Rabbit,IgG                                                                                                                                                                                                                                                                                                                                       |
| <b>Purification</b>                | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                                                                        |
| <b>Dilution</b>                    | WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/10000.. IF 1:50-200                                                                                                                                                                                                                                                                                         |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                      |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                                                                                                                                                                         |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                 |
| <b>Synonyms</b>                    | DBP5;DDX19;DDX19B;DEAD box protein 19B;DEAD box RNA helicase DEAD5; RNAh;ATP-dependent RNA helicase DDX19B;TDBP                                                                                                                                                                                                                                              |
| <b>Observed Band</b>               | 50kD                                                                                                                                                                                                                                                                                                                                                         |
| <b>Calculated Molecular Weight</b> | 50kD                                                                                                                                                                                                                                                                                                                                                         |
| <b>Cell Pathway</b>                | Cytoplasm . Nucleus, nucleoplasm . Associates with the nuclear pore complex cytoplasmic fibrils. .                                                                                                                                                                                                                                                           |
| <b>Tissue Specificity</b>          | Brain,Liver,Ovary,Skin,Stomach,Testis,                                                                                                                                                                                                                                                                                                                       |
| <b>Function</b>                    | function:ATP-dependent RNA helicase involved in mRNA export from the nucleus.,similarity:Belongs to the DEAD box helicase family. DDX19/DBP5 subfamily.,similarity:Contains 1 helicase ATP-binding domain., similarity:Contains 1 helicase C-terminal domain.,subcellular location:Nuclear pore complex cytoplasmic fibrils.,subunit:Interacts with NUP214., |
| <b>Background</b>                  | DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome                                           |



assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This gene encodes a DEAD box protein, which exhibits RNA-dependent ATPase and ATP-dependent RNA-unwinding activities. This protein is recruited to the cytoplasmic fibrils of the nuclear pore complex, where it participates in the export of mRNA from the nucleus. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],

**matters needing attention**


Avoid repeated freezing and thawing!

**Usage suggestions**


This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images





Western Blot analysis of various cells using DDX19B Polyclonal Antibody



Western blot analysis of lysates from HT-29 cells, using DDX19B Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using DDX19B antibody.



Immunohistochemical analysis of paraffin-embedded human brain tumor. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).