




TMEM101 Polyclonal Antibody


Catalog No	YP-Ab-13703
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	TMEM101
Immunogen	The antiserum was produced against synthesized peptide derived from human TMEM101. AA range:101-150
Specificity	TMEM101 Polyclonal Antibody detects endogenous levels of TMEM101 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TMEM101;Transmembrane protein 101;Putative NF-kappa-B-activating protein 130
Observed Band	28kD
Calculated Molecular Weight	28kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Lung,Muscle,
Function	function:May activate NF-kappa-B signaling pathways.,
Background	function:May activate NF-kappa-B signaling pathways.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of TMEM101 Antibody. The lane on the right is blocked with the TMEM101 peptide.



Western blot analysis of the lysates from HepG2 cells using TMEM101 antibody.