



# CD168 Polyclonal Antibody

|                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-Ab-13761                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Reactivity</b>                  | Human;Rat;Mouse;                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Applications</b>                | WB;IHC;IF;ELISA                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Gene Name</b>                   | HMMR IHABP RHAMM                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Immunogen</b>                   | Synthetic peptide from human protein at AA range: 280-340                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Specificity</b>                 | The antibody detects endogenous CD168                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                                                                                                                                                                            |
| <b>Source</b>                      | Polyclonal, Rabbit,IgG                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Purification</b>                | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                                                                                                                                              |
| <b>Dilution</b>                    | WB 1:500-2000,IHC-p 1:500-200, ELISA 1:10000-20000. IF 1:50-200                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Synonyms</b>                    | CD168;IHABP;RHAMM;CD168/RHAMM;HMMR IHABP RHAMM;Hyaluronan mediated motility receptor (Intracellular hyaluronic acid-binding protein) (Receptor for hyaluronan-mediated motility) (CD antigen CD168);Hyaluronan mediated motility receptor (Intracellular hyaluronic acid-binding protein;Receptor for hyaluronan-mediated motility;CD antigen CD168)                                                                               |
| <b>Observed Band</b>               | 84kD                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Calculated Molecular Weight</b> | 84kD                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Cell Pathway</b>                | Cell surface . Cytoplasm . Cytoplasm, cytoskeleton, spindle .                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Tissue Specificity</b>          | Expressed in testis (PubMed:22965910). Expressed in the breast (PubMed:8890751).                                                                                                                                                                                                                                                                                                                                                   |
| <b>Function</b>                    | function:Involved in cell motility. When hyaluronan binds to HMMR, the phosphorylation of a number of proteins, including the focal adhesion kinase occurs. May also be involved in cellular transformation and metastasis formation, and in regulating extracellular-regulated kinase (ERK) activity.,subunit:Subunit of the HARC complex.,tissue specificity:Expressed in breast cancer cell lines and in normal breast tissue., |



## Background

The protein encoded by this gene is involved in cell motility. It is expressed in breast tissue and together with other proteins, it forms a complex with BRCA1 and BRCA2, thus is potentially associated with higher risk of breast cancer. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Dec 2008],


## matters needing attention

Avoid repeated freezing and thawing!


## Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.


## Products Images



Western blot analysis of 293T lysate, antibody was diluted at 1000. Secondary antibody(catalog#:RS0002) was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200