




A1BG Monoclonal Antibody


Catalog No	YP-Ab-13783
Isotype	IgG
Reactivity	Human
Applications	WB;IHC;IF;ELISA
Gene Name	A1BG
Immunogen	Purified recombinant fragment of human A1BG expressed in E. Coli.
Specificity	A1BG Monoclonal Antibody detects endogenous levels of A1BG protein.
Formulation	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	WB: 1/500 - 1/2000. IHC: 1/200 - 1/1000. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A1BG;Alpha 1B-Glycoprotein;A1B;ABG;alpha 1 B glycoprotein;Alpha-1B-glycoprotein;Alpha-1-B glycoprotein
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted.
Tissue Specificity	Plasma.
Function	online information:The Singapore human mutation and polymorphism database, similarity:Contains 5 Ig-like V-type (immunoglobulin-like) domains.,tissue specificity:Plasma.,
Background	The protein encoded by this gene is a plasma glycoprotein of unknown function. The protein shows sequence similarity to the variable regions of some immunoglobulin supergene family member proteins. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis using A1BG Monoclonal Antibody against A1BG-hlgGfc transfected HEK293 (3)cell lysate.



Immunohistochemistry analysis of paraffin-embedded human Kidney tissues with AEC staining using A1BG Monoclonal Antibody.