



Helios Polyclonal Antibody

Catalog No	YP-Ab-01770
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	IKZF2
Immunogen	The antiserum was produced against synthesized peptide derived from human IKZF2. AA range:331-380
Specificity	Helios Polyclonal Antibody detects endogenous levels of Helios protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IKZF2;ZNFN1A2;ZNF1A2;Zinc finger protein Helios;IKZF 2;Helios;Ikaros family zinc finger protein 2
Observed Band	56kD
Calculated Molecular Weight	56kD
Cell Pathway	Nucleus .
Tissue Specificity	Epithelium,Testis,
Function	function:Associates with Ikaros at centromeric heterochromatin., similarity:Belongs to the Ikaros C2H2-type zinc-finger protein family., similarity:Contains 6 C2H2-type zinc fingers.,subunit:Interacts with IKZF4 AND IKZF5.,
Background	This gene encodes a member of the Ikaros family of zinc-finger proteins. Three members of this protein family (Ikaros, Aiolos and Helios) are hematopoietic-specific transcription factors involved in the regulation of lymphocyte development. This protein forms homo- or hetero-dimers with other Ikaros family members, and is thought to function predominantly in early hematopoietic



development. Multiple transcript variants encoding different isoforms have been found for this gene, but the biological validity of some variants has not been determined. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of HeLa cells using Helios Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA).



Western blot analysis of the lysates from HeLa cells using IKZF2 antibody.