



NPAS4 Polyclonal Antibody

Catalog No	YP-Ab-01913
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	IHC;IF;ELISA
Gene Name	NPAS4
Immunogen	The antiserum was produced against synthesized peptide derived from human NPAS4. AA range:603-652
Specificity	NPAS4 Polyclonal Antibody detects endogenous levels of NPAS4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/20000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NPAS4;Neuronal PAS domain-containing protein 4;BHLHE79;NXF;PASD10; Neuronal PAS4;Class E basic helix-loop-helix protein 79;HLH-PAS transcription factor NXF;PAS domain-containing protein 10
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Brain.
Function	function:Acts as a transcriptional activator in the presence of ARNT. Can activate the CME (CNS midline enhancer) element and the expression of the drebrin gene.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain., similarity:Contains 1 PAC (PAS-associated C-terminal) domain., similarity:Contains 2 PAS (PER-ARNT-SIM) domains.,subunit:Efficient DNA binding requires dimerization with another bHLH protein. Forms a heterodimer with ARNT.,tissue specificity:Brain.,



Background

NXF is a member of the basic helix-loop-helix-PER (MIM 602260)-ARNT (MIM 126110)-SIM (see SIM2; MIM 600892) (bHLH-PAS) class of transcriptional regulators, which are involved in a wide range of physiologic and developmental events (Ooe et al., 2004 [PubMed 14701734]).[supplied by OMIM, Mar 2008],



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using NPAS4 Antibody. The picture on the right is blocked with the synthesized peptide.