




NRIP2 Polyclonal Antibody

Catalog No	YP-Ab-01918
Isotype	IgG
Reactivity	Human;Mouse
Applications	IHC;IF;ELISA
Gene Name	NRIP2
Immunogen	The antiserum was produced against synthesized peptide derived from human NRIP2. AA range:1-50
Specificity	NRIP2 Polyclonal Antibody detects endogenous levels of NRIP2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/20000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NRIP2;Nuclear receptor-interacting protein 2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Amygdala,Brain,PCR rescued clones,Trachea,
Function	function:Down-regulates transcriptional activation by nuclear receptors such as NR1F2.,subunit:Interacts with NR1F2, RARA and THRB in a ligand-dependent manner.,
Background	function:Down-regulates transcriptional activation by nuclear receptors such as NR1F2.,subunit:Interacts with NR1F2, RARA and THRB in a ligand-dependent manner.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using NRIP2 Antibody. The picture on the right is blocked with the synthesized peptide.