




NY-REN-50 Polyclonal Antibody

Catalog No	YP-Ab-01926
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	ZNF839
Immunogen	Synthesized peptide derived from the Internal region of human NY-REN-50.
Specificity	NY-REN-50 Polyclonal Antibody detects endogenous levels of NY-REN-50 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Zinc finger protein 839;Renal carcinoma antigen NY-REN-50;C14orf131;NY-REN-50;ZNF839
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Amygdala,Lung,Neuroblastoma,Placenta,PNS,Testis,Thymus,
Function	similarity:Contains 1 C2H2-type zinc finger.,
Background	similarity:Contains 1 C2H2-type zinc finger.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).