




RECS1 Polyclonal Antibody


Catalog No	YP-Ab-01972
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;IHC;IF;ELISA
Gene Name	TMBIM1
Immunogen	The antiserum was produced against synthesized peptide derived from human TMBIM1. AA range:251-300
Specificity	RECS1 Polyclonal Antibody detects endogenous levels of RECS1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RECS1;TMBIM1;Protein lifeguard 3;LFG3;PP1201;PSEC0158;Protein RECS1 homolog;Transmembrane BAX inhibitor motif-containing protein 1
Observed Band	35kD
Calculated Molecular Weight	35kD
Cell Pathway	Membrane ; Multi-pass membrane protein . Lysosome membrane . Endosome membrane .
Tissue Specificity	Aorta,Brain,Colon,Pancreas,Placenta,Platelet,Tongue,Trachea
Function	similarity:Belongs to the BI1 family.,
Background	similarity:Belongs to the BI1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images




Western Blot analysis of various cells using RECS1 Polyclonal Antibody diluted at 1:1000



Western blot analysis of lysates from HUVEC cells, using TMBIM1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from K562 cells using TMBIM1 antibody.



Immunohistochemical analysis of paraffin-embedded human colon cancer. 1, Tris-EDTA, pH9.0 was used for antigen retrieval. 2 Antibody was diluted at 1:200 (4° overnight). 3, Secondary antibody was diluted at 1:200 (room temperature, 45min).