



RNF149 Polyclonal Antibody

Catalog No	YP-Ab-01988
Isotype	IgG
Reactivity	Human;Mouse
Applications	IHC;IF;ELISA
Gene Name	RNF149
Immunogen	The antiserum was produced against synthesized peptide derived from human RNF149. AA range:211-260
Specificity	RNF149 Polyclonal Antibody detects endogenous levels of RNF149 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RNF149;E3 ubiquitin-protein ligase RNF149;DNAPTP2;DNA polymerase-transactivated protein 2;RING finger protein 149
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	Brain,Cervix,Placenta,Teratocarcinoma,
Function	domain:The RING-type zinc finger domain mediates binding to an E2 ubiquitin-conjugating enzyme.,pathway:Protein modification; protein ubiquitination.,similarity:Contains 1 PA (protease associated) domain.,similarity:Contains 1 RING-type zinc finger.,
Background	domain:The RING-type zinc finger domain mediates binding to an E2 ubiquitin-conjugating enzyme.,pathway:Protein modification; protein ubiquitination.,similarity:Contains 1 PA (protease associated) domain.,similarity:Contains 1 RING-type zinc finger.,



**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using RNF149 Antibody. The picture on the right is blocked with the synthesized peptide.