



# UBF Polyclonal Antibody

|                                    |                                                                                                                                                                                                                                                                                                                                                                                         |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-Ab-02147                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Reactivity</b>                  | Human;Mouse;Rat                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Applications</b>                | WB;ELISA                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Gene Name</b>                   | UBTF                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Immunogen</b>                   | Synthesized peptide derived from UBF . at AA range: 420-500                                                                                                                                                                                                                                                                                                                             |
| <b>Specificity</b>                 | UBF Polyclonal Antibody detects endogenous levels of UBF protein.                                                                                                                                                                                                                                                                                                                       |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                                                                                                                                 |
| <b>Source</b>                      | Polyclonal, Rabbit,IgG                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Purification</b>                | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                                                                                                   |
| <b>Dilution</b>                    | Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.                                                                                                                                                                                                                                                                                                     |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Synonyms</b>                    | UBTF;UBF-1;UBF 1;UBF;NOR 90;Nucleolar transcription factor 1;UBF1; Autoantigen NOR-90;Upstream-binding factor 1                                                                                                                                                                                                                                                                         |
| <b>Observed Band</b>               | 85-90 kDa                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Calculated Molecular Weight</b> | 89 kDa                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cell Pathway</b>                | Nucleus, nucleolus .                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Tissue Specificity</b>          | Lymph,Placenta,Testis,                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Function</b>                    | function:Recognizes the ribosomal RNA gene promoter and activates transcription mediated by RNA polymerase I through cooperative interactions with the species-specific factor SL1. It binds specifically to the upstream control element.,PTM:Phosphorylated and activated by PIK3CA.,similarity:Contains 6 HMG box DNA-binding domains.,subunit:Homodimer. Binds to IRS1 and PIK3CA., |
| <b>Background</b>                  | This gene encodes a member of the HMG-box DNA-binding protein family. The encoded protein plays a critical role in ribosomal RNA transcription as a key component of the pre-initiation complex, mediating the recruitment of RNA polymerase I to rDNA promoter regions. The encoded protein may also play                                                                              |



important roles in chromatin remodeling and pre-rRNA processing, and its activity is regulated by both phosphorylation and acetylation. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. Pseudogenes of this gene are located on the short arm of chromosomes 3, 11 and X and the long arm of chromosome 11. [provided by RefSeq, Aug 2011],

**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images