




ABHD12B Polyclonal Antibody

Catalog No	YP-Ab-02452
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	ABHD12B
Immunogen	The antiserum was produced against synthesized peptide derived from human ABHD12B. AA range:281-330
Specificity	ABHD12B Polyclonal Antibody detects endogenous levels of ABHD12B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ABHD12B;BEM46L3;c14_5314;C14orf29;Abhydrolase domain-containing protein 12B
Observed Band	40kD
Calculated Molecular Weight	40kD
Cell Pathway	
Tissue Specificity	Kidney,Skin,
Function	similarity:Belongs to the serine esterase family.,
Background	similarity:Belongs to the serine esterase family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma, using ABHD12B Antibody. The picture on the right is blocked with the synthesized peptide.