



NOS2 Polyclonal Antibody

Catalog No	YP-Ab-02719
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	IF/ICC;WB;IHC-p;ELISA
Gene Name	NOS2, INOS
Immunogen	The antiserum was produced against synthesized peptide derived from human iNOS. AA range:117-166
Specificity	NOS2 Polyclonal Antibody detects endogenous levels of NOS2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IF: 1:50-200 WB 1:500-2000, ELISA 1:10000-20000 IHC 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Molecular Weight	13KD
Observed Band	13KD
Cell Pathway	Cytoplasm, cytosol . Localizes as discrete foci scattered throughout the cytosol and in the presence of SPSB1 and SPSB4, exhibits a more diffuse cytosolic localization.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images