



RN-tre Polyclonal Antibody

Catalog No	YP-Ab-02774
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	USP6NL
Immunogen	The antiserum was produced against synthesized peptide derived from human USP6NL. AA range:451-500
Specificity	RN-tre Polyclonal Antibody detects endogenous levels of RN-tre protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	USP6 N-terminal-like protein;Related to the N-terminus of tre;KIAA0019;RN-tre; USP6NL
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus. Cytoplasmic vesicle.
Tissue Specificity	Widely expressed.
Function	function:Acts as a GTPase-activating protein for RAB5A. Involved in receptor trafficking. In complex with EPS8 inhibits internalization of EGFR., PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 Rab-GAP TBC domain.,subunit:Interacts with EPS8.,tissue specificity:Widely expressed.,
Background	function:Acts as a GTPase-activating protein for RAB5A. Involved in receptor trafficking. In complex with EPS8 inhibits internalization of EGFR., PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 Rab-GAP TBC domain.,subunit:Interacts with EPS8.,tissue specificity:Widely expressed.,


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using USP6NL Antibody. The picture on the right is blocked with the synthesized peptide.