





UBTD1 Polyclonal Antibody

Catalog No	YP-Ab-02820
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	IHC;IF;ELISA
Gene Name	UBTD1
Immunogen	The antiserum was produced against synthesized peptide derived from human UBTD1. AA range:151-200
Specificity	UBTD1 Polyclonal Antibody detects endogenous levels of UBTD1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/20000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ubiquitin domain containing 1;UBTD1;Ubiquitin domain-containing protein 1
Observed Band	26kD
Calculated Molecular Weight	26kD
Cell Pathway	
Tissue Specificity	Coronary artery,Embryo,Placenta,
Function	similarity:Contains 1 ubiquitin-like domain.,
Background	The degradation of many proteins is carried out by the ubiquitin pathway in which proteins are targeted for degradation by covalent conjugation of the polypeptide ubiquitin. This gene encodes a protein that belongs to the ubiquitin family of proteins. The encoded protein is thought to regulate E2 ubiquitin conjugating enzymes belonging to the UBE2D family. [provided by RefSeq, Mar 2014],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using UBTD1 Antibody. The picture on the right is blocked with the synthesized peptide.