




USP32 Polyclonal Antibody

Catalog No	YP-Ab-02833
Isotype	IgG
Reactivity	Human;Mouse
Applications	IHC;IF;WB;ELISA
Gene Name	USP32
Immunogen	The antiserum was produced against synthesized peptide derived from human USP32. AA range:271-320
Specificity	USP32 Polyclonal Antibody detects endogenous levels of USP32 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	USP10;Ubiquitin-specific-processing protease 32;Ubiquitin carboxyl-terminal hydrolase 32;Renal carcinoma antigen NY-REN-60;NY REN 60;USP32;Deubiquitinating enzyme 32;Ubiquitin thioesterase 32
Observed Band	180kD
Calculated Molecular Weight	180kD
Cell Pathway	Membrane ; Lipid-anchor . Golgi apparatus .
Tissue Specificity	Brain,Epithelium,Peripheral Nervous System,Skin,Testis,Uterus,
Function	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., similarity:Belongs to the peptidase C19 family.,similarity:Contains 1 DUSP domain.,similarity:Contains 3 EF-hand domains.,
Background	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., similarity:Belongs to the peptidase C19 family.,similarity:Contains 1 DUSP domain.,similarity:Contains 3 EF-hand domains.,
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using USP32 Antibody. The picture on the right is blocked with the synthesized peptide.