



# UHRF1 Polyclonal Antibody

|                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-Ab-02848                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Reactivity</b>                  | Human;Rat;Mouse;                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Applications</b>                | WB;ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Gene Name</b>                   | UHRF1                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Immunogen</b>                   | Synthesized peptide derived from the Internal region of human UHRF1.                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Specificity</b>                 | UHRF1 Polyclonal Antibody detects endogenous levels of UHRF1 protein.                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Source</b>                      | Polyclonal, Rabbit,IgG                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Purification</b>                | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                                                                                                                                                                                          |
| <b>Dilution</b>                    | Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Synonyms</b>                    | ICBP90;HuNp95;hUHRF1;hNp95;EC:2.3.2.27;UHRF1;E3 ubiquitin-protein ligase UHRF1;NP95;RNF106;Inverted CCAAT box-binding protein of 90 kDa; Nuclear protein 95;Nuclear zinc finger protein Np95;RING finger protein 106; Transcription factor ICBP90;Ubiquitin-like PHD and RING finger domain-containing protein 1;Ubiquitin-like-containing PHD and RING finger domains protein 1                                                                                               |
| <b>Observed Band</b>               | 89kD                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Calculated Molecular Weight</b> | 89kD                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Cell Pathway</b>                | Nucleus . Localizes to replication foci. Enriched in pericentric heterochromatin. Also localizes to euchromatic regions.                                                                                                                                                                                                                                                                                                                                                       |
| <b>Tissue Specificity</b>          | Expressed in thymus, bone marrow, testis, lung and heart. Overexpressed in breast cancer.                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Function</b>                    | developmental stage:Expressed in fetal thymus, liver and kidney.,domain:The RING finger is required for ubiquitin ligase activity.,domain:The YDG domain mediates the interaction with histone H3.,function:Putative E3 ubiquitin-protein ligase. May participate in methylation-dependent transcriptional regulation. Binds to inverted 5'-CCAAT-3' box 2 in the TOP2A promoter, and activates TOP2A expression. Important for G1/S transition. May be involved in DNA repair |



and chromosomal stability.,induction:Up-regulated in proliferating cells, and down-regulated in quiescent cells. Down-regulated upon adriamycin-induced DNA damage, in a TP53/p53 and CDKN1A-dependent way. Induced by E2F1 transcription factor.,pathway:Protein modification; protein ubiquitination., PTM:Phosphorylated on serine residues. Phosphorylation may enhance DNA-binding activity.,PTM:Ubiquitinated; which leads to proteasomal degrad

## Background

This gene encodes a member of a subfamily of RING-finger type E3 ubiquitin ligases. The protein binds to specific DNA sequences, and recruits a histone deacetylase to regulate gene expression. Its expression peaks at late G1 phase and continues during G2 and M phases of the cell cycle. It plays a major role in the G1/S transition by regulating topoisomerase IIalpha and retinoblastoma gene expression, and functions in the p53-dependent DNA damage checkpoint. It is regarded as a hub protein for the integration of epigenetic information. This gene is up-regulated in various cancers, and it is therefore considered to be a therapeutic target. Multiple transcript variants encoding different isoforms have been found for this gene. A related pseudogene exists on chromosome 12. [provided by RefSeq, Feb 2014],


## matters needing attention

Avoid repeated freezing and thawing!

## Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images



Western Blot analysis of extracts from Jurkat cells, using UHRF1 Polyclonal Antibody. Secondary antibody(catalog#:RS0002) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).