



MYL2 Monoclonal Antibody

Catalog No	YP-Ab-02936
Isotype	IgG
Reactivity	Human
Applications	WB;ELISA
Gene Name	MYL2
Immunogen	Purified recombinant fragment of MYL2 expressed in E. Coli.
Specificity	MYL2 Monoclonal Antibody detects endogenous levels of MYL2 protein.
Formulation	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MLC 2v;MYL2;MLC 2;MLC2;MLC-2;Myosin regulatory light chain 2 ventricular/cardiac muscle isoform;Myosin regulatory light chain 2;ventricular/cardiac muscle isoform;MLC-2v
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, myofibril, sarcomere, A band .
Tissue Specificity	Highly expressed in type I muscle fibers.
Function	disease:Defects in MYL2 are the cause of cardiomyopathy familial hypertrophic type 10 (CMH10) [MIM:608758]. Familial hypertrophic cardiomyopathy is a hereditary heart disorder characterized by ventricular hypertrophy, which is usually asymmetric and often involves the interventricular septum. The symptoms include dyspnea, syncope, collapse, palpitations, and chest pain. They can be readily provoked by exercise. The disorder has inter- and intrafamilial variability ranging from benign to malignant forms with high risk of cardiac failure and sudden cardiac death.,disease:Defects in MYL2 are the cause of cardiomyopathy hypertrophic with mid-left ventricular chamber type 2 (MVC2) [MIM:608758]. MVC2 is a very rare variant of familial hypertrophic cardiomyopathy, characterized by mid-left ventricular chamber thickening.,



miscellaneous: This chain binds calcium., similarity: Contains 3 EF-hand domain

Background

Thus gene encodes the regulatory light chain associated with cardiac myosin beta (or slow) heavy chain. Ca²⁺ triggers the phosphorylation of regulatory light chain that in turn triggers contraction. Mutations in this gene are associated with mid-left ventricular chamber type hypertrophic cardiomyopathy. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis using MYL2 Monoclonal Antibody against rat fetal heart tissue lysate.