



MYL3 Monoclonal Antibody

Catalog No	YP-Ab-02937
Isotype	IgG
Reactivity	Human
Applications	WB;IHC;IF;ELISA
Gene Name	MYL3
Immunogen	Purified recombinant fragment of MYL3 expressed in E. Coli.
Specificity	MYL3 Monoclonal Antibody detects endogenous levels of MYL3 protein.
Formulation	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	WB: 1/500 - 1/2000. IHC: 1/200 - 1/1000. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MLC1V;MLC-IV/sb;MYL 3;Myosin light chain 1;slow-twitch muscle B/ventricular isoform;Myosin light chain 3;MYL3;Cardiac myosin light chain 1;CMLC1;MLC1SB;Ventricular/slow twitch myosin alkali light chain
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytosol,muscle myosin complex,myosin complex,sarcomere,A band,I band,
Tissue Specificity	Heart,Skeletal muscle,
Function	disease:Defects in MYL3 are the cause of cardiomyopathy familial hypertrophic type 8 (CMH8) [MIM:608751]. Familial hypertrophic cardiomyopathy is a hereditary heart disorder characterized by ventricular hypertrophy, which is usually asymmetric and often involves the interventricular septum. The symptoms include dyspnea, syncope, collapse, palpitations, and chest pain. They can be readily provoked by exercise. The disorder has inter- and intrafamilial variability ranging from benign to malignant forms with high risk of cardiac failure and sudden cardiac death. CMH8 inheritance can be autosomal dominant or recessive.,disease:Defects in MYL3 are the cause of cardiomyopathy hypertrophic with mid-left ventricular chamber type 1 (MVC1) [MIM:608751]. MVC1 is a very rare variant of familial hypertrophic



cardiomyopathy, characterized by mid-left ventricular chamber thickening.,
function:Regulatory

Background

MYL3 encodes myosin light chain 3, an alkali light chain also referred to in the literature as both the ventricular isoform and the slow skeletal muscle isoform. Mutations in MYL3 have been identified as a cause of mid-left ventricular chamber type hypertrophic cardiomyopathy. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis using MYL3 Monoclonal Antibody against rat fetal heart tissues lysate.



Immunohistochemistry analysis of paraffin-embedded human Heart tissues with AEC staining using MYL3 Monoclonal Antibody.