



β I tubulin Monoclonal Antibody

Catalog No	YP-Ab-02984
Isotype	IgG
Reactivity	Zebrafish
Applications	WB
Gene Name	TUBB1
Immunogen	Synthetic Peptide of β I tubulin
Specificity	The antibody detects Zebrafish endogenous β I tubulin proteins.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	WB: 1:500-10000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Tubulin β;TUBB1;Tubulin beta-1 chain;TUBB2A;TUBB2;Tubulin beta-2A chain;TUBB2B;Tubulin beta-2B chain;TUBB4B;TUBB2C;Tubulin beta-4B chain;Tubulin beta-2 chain;Tubulin beta-2C chain;TUBB3;TUBB4;Tubulin beta-3 chain;Tubulin beta-4 chain;Tubulin beta-III;TUBB4A;TUBB5;Tubulin beta-4A chain;Tubulin 5 beta;TUBB;OK/SW-cl.56;Tubulin beta chain;Tubulin beta-5 chain;TUBB6;Tubulin beta-6 chain;TUBB8;Tubulin beta-8 chain
Observed Band	50kD
Calculated Molecular Weight	50kD
Cell Pathway	Cytoplasm, cytoskeleton .
Tissue Specificity	Hematopoietic cell-specific. Major isotype in leukocytes, where it represents 50% of all beta-tubulins.
Function	function:Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha-chain.,similarity:Belongs to the tubulin family.,subunit:Dimer of alpha and beta chains.,



Background

This gene encodes a member of the beta tubulin protein family. Beta tubulins are one of two core protein families (alpha and beta tubulins) that heterodimerize and assemble to form microtubules. This protein is specifically expressed in platelets and megakaryocytes and may be involved in proplatelet production and platelet release. A mutations in this gene is associated with autosomal dominant macrothrombocytopenia. Two pseudogenes of this gene are found on chromosome Y.[provided by RefSeq, Jul 2010],

matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

94KD
66KD
45KD
35KD
26KD
20KD



Western blot analysis of Zebrafish skeletal muscle, diluted at 1:5000.