



POTE15 Polyclonal Antibody

Catalog No	YP-Ab-03172
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;IHC;IF;ELISA
Gene Name	POTEB
Immunogen	The antiserum was produced against synthesized peptide derived from human A26B1. AA range:481-530
Specificity	POTE15 Polyclonal Antibody detects endogenous levels of POTE15 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/5000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	POTEB;A26B1;ANKRD26 like family B member 1;POTE 15;POTE ankyrin domain family member B3;POTE15;POTE ankyrin domain family member B; ANKRD26-like family B member 1;Prostate;ovary;testis-expressed protein on chromosome 15;POTE-15
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	
Tissue Specificity	
Function	
Background	POTEB (POTE Ankyrin Domain Family Member B) is a Protein Coding gene. An important paralog of this gene is ANKRD18A.
matters needing attention	Avoid repeated freezing and thawing!




Usage suggestions


This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images


Mouse-kidney 3T3



Western blot analysis of Mouse-kidney 3T3 lysis using POTE15 antibody. Antibody was diluted at 1:500



Western blot analysis of lysates from Jurkat cells, using A26B1 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).