




Nuclear Receptor NR4A1 (Phospho-Ser351) Antibody

Catalog No	YP-Ab-03301
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB; ELISA
Gene Name	NR4A1 GFRP1 HMR NAK1
Immunogen	Synthesized pospho peptide derived from human Nuclear Receptor NR4A1 (Phospho-Ser351)
Specificity	This antibody detects endogenous pospho levels of human Nuclear Receptor NR4A1 (Phospho-Ser351)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000, ELISA(peptide)1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Early response protein NAK1;GFRP 1;GFRP;GFRP1;Growth factor inducible nuclear protein N1;Growth Factor Inducible Nuclear Protein NP1;Growth Factor Response Protein 1;Hbr1;HMR;Hormone Receptor;MGC9485;N1;N1 nuclear protein;NAK 1;NAK1;Nerve growth factor IB nuclear receptor variant 1;NGFIB; NP 1;NP1;NR4A1;NR4A1_HUMAN;Nuclear hormone receptor NUR/77;Nuclear Hormone Receptor TR3;Nuclear receptor subfamily 4 group A member 1; NUR77;mouse;homolog of;Orphan nuclear receptor HMR;Orphan nuclear receptor NR4A1;Orphan nuclear receptor TR3;Orphan receptor tr3;Receptor NGFIB;ST 59;ST-59;ST59;Steroid receptor TR3;Testicular receptor 3;TR 3;TR3; TR3 orphan receptor;NR4A1 GFRP1 HMR NAK1;Nuclear Receptor NR4A1 (Phospho-Ser351);Nuclear receptor subfamily 4 group A member 1 (Early response protein NAK1;Testicular receptor 3)
Observed Band	77kD
Calculated Molecular Weight	77kD
Cell Pathway	Nucleus . Cytoplasm . Mitochondrion . Nuclear export to the cytoplasm is XPO1-mediated and positively regulated by IFI27 (PubMed:22427340). Translocation to the mitochondrion upon interaction with RXRA and upon the presence of 9-cis retinoic acid (PubMed:17761950). .



Tissue Specificity	Fetal muscle and adult liver, brain and thyroid.
Function	function:Orphan nuclear receptor.,induction:By growth-stimulating agents.,similarity:Belongs to the nuclear hormone receptor family.,similarity:Belongs to the nuclear hormone receptor family. NR4 subfamily.,similarity:Contains 1 nuclear receptor DNA-binding domain.,subunit:Interacts with GADD45GIP1.,tissue specificity:Fetal muscle and adult liver, brain and thyroid.,
Background	This gene encodes a member of the steroid-thyroid hormone-retinoid receptor superfamily. Expression is induced by phytohemagglutinin in human lymphocytes and by serum stimulation of arrested fibroblasts. The encoded protein acts as a nuclear transcription factor. Translocation of the protein from the nucleus to mitochondria induces apoptosis. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2011],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of various lysates, primary antibody was diluted at 1:1000, 4° over night, secondary antibody(cat: RS23920)was diluted at 1:10000, 37° 1hour.