



RBBP5 mouse mAb

Catalog No	YP-Ab-03477
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;IF
Gene Name	rbbp5
Immunogen	Purified recombinant human RBBP5 protein fragments expressed in E.coli.
Specificity	This antibody detects endogenous levels of RBBP5 and does not cross-react with related proteins.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution	wb dilution 1:1000 icc dilution 1:100. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RBQ3;SWD1;RBBP5;RBBP 5;RBBP-5;RBBP5_HUMAN;RBQ 3;Retinoblastoma binding protein 5;Retinoblastoma binding protein RBQ3;Retinoblastoma binding protein5;Retinoblastoma-binding protein 5;Retinoblastoma-binding protein RBQ-3;Set1c WD40 repeat protein;homolog.
Observed Band	70kD
Calculated Molecular Weight	70kD
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitously expressed.
Function	function: Binds preferentially to underphosphorylated retinoblastoma protein., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Contains 6 WD repeats., subunit: Component of the SET1 complex, at least composed of the catalytic subunit (SETD1A or SETD1B), WDR5, WDR82, RBBP5, ASH2/ASH2L and CXXC1/CFP1. Component of MLL-containing complexes (named MLL, ASCOM, MLL2/MLL3 or MLL3/MLL4 complex): at least composed ASH2L, RBBP5, DPY30, WDR5, one or several histone methyltransferases (MLL, MLL2, MLL3 and/or MLL4), and the facultative components MEN1, HCFC1, HCFC2, NCOA6, KDM6A, PAXIP1/PTIP and



C16orf53/PA1.,tissue specificity:Ubiquitously expressed.,

Background

This gene encodes a ubiquitously expressed nuclear protein which belongs to a highly conserved subfamily of WD-repeat proteins. The encoded protein binds directly to retinoblastoma protein, which regulates cell proliferation. It interacts preferentially with the underphosphorylated retinoblastoma protein via the E1A-binding pocket B. Three alternatively spliced transcript variants that encode different protein isoforms have been described for this gene. [provided by RefSeq, Jul 2010],


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot detection of RBBP5 in C2C12,THP-1,PC12 and 3T3 cell lysates using RBBP5 mouse mAb (1:1000 diluted).Predicted band size:70KDa.Observed band size:70KDa.



Immunofluorescent analysis of HeLa cells fixed with 4% Paraformaldehyde and using anti-RBBP5 mouse mAb (dilution 1:100). DAPI was used to stain nucleus(blue).