



ANKRD20 Polyclonal Antibody


Catalog No	YP-Ab-03698
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	ANKRD20A3
Immunogen	The antiserum was produced against synthesized peptide derived from human ANKRD20A3. AA range:551-600
Specificity	ANKRD20 Polyclonal Antibody detects endogenous levels of ANKRD20 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A2A3_HUMAN;ANKRD2A3;Ankyrin repeat domain 2 family member A3;Ankyrin repeat domain 2A related;Ankyrin repeat domain containing protein 2A3;Ankyrin repeat domain-containing protein 2A3;MGC176486;MGC19858;ANKRD20; ANKRD20A3;Ankyrin repeat domain-containing protein 20A3
Observed Band	95kD
Calculated Molecular Weight	95kD
Cell Pathway	
Tissue Specificity	Testis,
Function	similarity:Contains 5 ANK repeats.,tissue specificity:Widely expressed; not found in kidney, nor in liver.,
Background	similarity:Contains 5 ANK repeats.,tissue specificity:Widely expressed; not found in kidney, nor in liver.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemistry analysis of paraffin-embedded human brain, using ANKRD20A3 Antibody. The picture on the right is blocked with the synthesized peptide.