



Annexin X Polyclonal Antibody


Catalog No	YP-Ab-03705
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	ANXA10
Immunogen	The antiserum was produced against synthesized peptide derived from human ANXA10. AA range:261-310
Specificity	Annexin X Polyclonal Antibody detects endogenous levels of Annexin X protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Annexin A10;Annexin 10;Annexin 14;ANX10;ANX14;Annexin X;ANXA10; Annexin-10;Annexin-14
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	mitochondrion,
Tissue Specificity	Colon,Lung,
Function	similarity:Belongs to the annexin family.,similarity:Contains 4 annexin repeats.,
Background	This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The function of this gene has not yet been determined. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!




Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Annexin X Polyclonal Antibody



Western blot analysis of the lysates from Jurkat cells using ANXA10 antibody.