




eIF1B Polyclonal Antibody

Catalog No	YP-Ab-03836
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	IHC;IF;ELISA
Gene Name	EIF1B
Immunogen	Synthesized peptide derived from the N-terminal region of human eIF1B. AA range: 26-75
Specificity	eIF1B Polyclonal Antibody detects endogenous levels of eIF1B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/20000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	eIF1b;EIF1B_HUMAN;Eukaryotic translation initiation factor 1b;Protein translation factor SUI1 homolog GC2;Protein translation factor SUI1 homolog GC20
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Blood,Cerebellum,Cerebrum,Kidney,Liver,Pituitary,
Function	function:Probably involved in translation.,similarity:Belongs to the SUI1 family.,
Background	function:Probably involved in translation.,similarity:Belongs to the SUI1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).