




FAM3D Polyclonal Antibody

Catalog No	YP-Ab-03867
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;IHC;IF;ELISA
Gene Name	FAM3D
Immunogen	Synthesized peptide derived from the C-terminal region of human FAM3D.
Specificity	FAM3D Polyclonal Antibody detects endogenous levels of FAM3D protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UNQ567/PRO1130;Protein FAM3D;OIT1;EF7;FAM3D
Observed Band	25kD
Calculated Molecular Weight	25kD
Cell Pathway	Secreted .
Tissue Specificity	Abundantly expressed in placenta and weakly expressed in small intestine.
Function	similarity:Belongs to the FAM3 family.,tissue specificity:Abundantly expressed in placenta and weakly expressed in small intestine.,
Background	similarity:Belongs to the FAM3 family.,tissue specificity:Abundantly expressed in placenta and weakly expressed in small intestine.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight). 2, Tris-EDTA, pH 9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200 (room temperature, 30 min).