




Gab 2/4 Polyclonal Antibody

Catalog No	YP-Ab-03885
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;IHC;IF;ELISA
Gene Name	GAB2/GAB4
Immunogen	The antiserum was produced against synthesized peptide derived from human GAB2. AA range:21-70
Specificity	Gab 2/4 Polyclonal Antibody detects endogenous levels of Gab 2/4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Gab 2/4;GAB2/GAB4;GRB2-associated-binding protein 2/4;GAB4;GRB2-associated-binding protein 4;GRB2-associated binder 2-like;GAB2-like;GRB2-associated binder 4;GRB2-associated-binding protein 2-like;Growth factor receptor bound protein 2-associated protein 4;GAB2;KIAA0571;GRB2-associate
Observed Band	62kD
Calculated Molecular Weight	62kD
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the GAB family.,similarity:Contains 1 PH domain.,
Background	similarity:Belongs to the GAB family.,similarity:Contains 1 PH domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemistry analysis of paraffin-embedded human esophagus, using GAB4 Antibody. The picture on the right is blocked with the synthesized peptide.