




Galactosidase β 1L3 Polyclonal Antibody

Catalog No	YP-Ab-03887
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	GLB1L3
Immunogen	The antiserum was produced against synthesized peptide derived from human BGLL3. AA range:291-340
Specificity	Galactosidase β 1L3 Polyclonal Antibody detects endogenous levels of Galactosidase β 1L3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	Galactosidase β 1L3;GLB1L3;Beta-galactosidase-1-like protein 3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	vacuole,
Tissue Specificity	Lung,Testis,
Function	similarity:Belongs to the glycosyl hydrolase 35 family.,
Background	similarity:Belongs to the glycosyl hydrolase 35 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma, using GLB1L3 Antibody. The picture on the right is blocked with the synthesized peptide.