



GCNT3 Polyclonal Antibody

Catalog No	YP-Ab-03892
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	GCNT3
Immunogen	The antiserum was produced against synthesized peptide derived from human GCNT3. AA range:226-275
Specificity	GCNT3 Polyclonal Antibody detects endogenous levels of GCNT3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Beta-1;3-galactosyl-O-glycosyl-glycoprotein beta-1;6-N-acetylglucosaminyltransferase 3;C2/4GnT;C2GnT mucin type;C2GnT-M;C2GnT-mucin type;GCNT3;hC2GnT-M;Core 2/core 4 beta-1;6-N-acetylglucosaminyltransferase
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Single-pass type II membrane protein .
Tissue Specificity	Primarily expressed in mucus-secreting tissues. Expressed in colon, kidney, small intestine, trachea, and stomach, where mucin is produced.
Function	catalytic activity:UDP-N-acetyl-D-glucosamine + beta-D-galactosyl-1,3-N-acetyl-D-galactosaminyl-R = UDP + beta-D-galactosyl-1,3-(N-acetyl-beta-D-glucosaminyl-1,6)-N-acetyl-D-galactosaminyl-R.,catalytic activity:UDP-N-acetyl-D-glucosamine + beta-D-galactosyl-1,4-N-acetyl-D-glucosaminyl-R = UDP + N-acetyl-beta-D-glucosaminyl-1,6-beta-D-galactosyl-1,4-N-acetyl-D-glucosaminyl-R.,function:Glycosyltransferase that can synthesize all known mucin beta 6 N-acetylglucosaminides. Mediates core 2 and core 4 O-glycan branching, 2 important steps in mucin-type biosynthesis. Has also I-branching enzyme



activity by converting linear into branched poly-N-acetyllactosaminoglycans, leading to introduce the blood group I antigen during embryonic development., induction: By all-trans retinoic acid (ATRA), TNF-alpha and IL13. Strongly down-regulated in colorectal cancer., online information: Core 2/core 4 beta-1,

Background

This gene encodes a member of the N-acetylglucosaminyltransferase family. The encoded protein is a beta-6-N-acetylglucosamine-transferase that catalyzes the formation of core 2 and core 4 O-glycans on mucin-type glycoproteins.[provided by RefSeq, Apr 2009],


matters needing attention

Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunofluorescence analysis of A549 cells, using GCNT3 Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue, using GCNT3 Antibody. The picture on the right is blocked with the synthesized peptide.