



IREβ Polyclonal Antibody

Catalog No	YP-Ab-03932
Isotype	IgG
Reactivity	Human;Mouse
Applications	IHC;IF;ELISA
Gene Name	ERN2
Immunogen	The antiserum was produced against synthesized peptide derived from human ERN2. AA range:841-890
Specificity	IRE1α/β Polyclonal Antibody detects endogenous levels of IRE1α/β protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IREβ;ERN2;Serine/threonine-protein kinase/endoribonuclease IRE2;IRE2; Endoplasmic reticulum-to-nucleus signaling 2;Inositol-requiring protein 2;hIRE2p; Ire1-beta;IRE1b
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type I membrane protein .
Tissue Specificity	Fetal brain,Fetal lung fibroblast,
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein., cofactor:Magnesium.,enzyme regulation:The kinase domain is activated by trans-autophosphorylation. Kinase activity is required for activation of the endoribonuclease domain.,function:Induces translational repression through 28S ribosomal RNA cleavage in response to ER stress. Pro-apoptotic. Appears to play no role in the unfolded-protein response, unlike closely related proteins., PTM:Autophosphorylated.,similarity:Belongs to the protein kinase superfamily. Ser/Thr protein kinase family.,similarity:Contains 1 KEN domain., similarity:Contains 1 protein kinase domain.,



Background

catalytic activity: ATP + a protein = ADP + a phosphoprotein.,
cofactor: Magnesium., enzyme regulation: The kinase domain is activated by trans-autophosphorylation. Kinase activity is required for activation of the endoribonuclease domain., function: Induces translational repression through 28S ribosomal RNA cleavage in response to ER stress. Pro-apoptotic. Appears to play no role in the unfolded-protein response, unlike closely related proteins., PTM: Autophosphorylated., similarity: Belongs to the protein kinase superfamily. Ser/Thr protein kinase family., similarity: Contains 1 KEN domain., similarity: Contains 1 protein kinase domain.,


matters needing attention

Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunofluorescence analysis of A549 cells, using ERN2 Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using ERN2 Antibody. The picture on the right is blocked with the synthesized peptide.