



KLHL29 Polyclonal Antibody

Catalog No	YP-Ab-03943
Isotype	IgG
Reactivity	Human;Mouse
Applications	IHC;IF;ELISA
Gene Name	KLHL29
Immunogen	The antiserum was produced against synthesized peptide derived from human KLHL29. AA range:261-310
Specificity	KLHL29 Polyclonal Antibody detects endogenous levels of KLHL29 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KBTBD9;Kelch like protein 29;KIAA1921;KLHL29;Kelch-like protein 29;Kelch repeat and BTB domain-containing protein 9
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cul3-RING ubiquitin ligase complex,
Tissue Specificity	Brain,Melanoma,Placenta,Spleen,
Function	caution:Although the complete sequence is not known with certainty, sequence shown here appears to be the most probable in accordance with the mouse sequence ortholog.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 6 Kelch repeats.,
Background	caution:Although the complete sequence is not known with certainty, sequence shown here appears to be the most probable in accordance with the mouse sequence ortholog.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 6 Kelch repeats.,



**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using KLHL29 Antibody. The picture on the right is blocked with the synthesized peptide.