



Ribosomal Protein L39L Polyclonal Antibody

Catalog No	YP-Ab-04147
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	RPL39L
Immunogen	The antiserum was produced against synthesized peptide derived from human RPL39L. AA range:1-50
Specificity	Ribosomal Protein L39L Polyclonal Antibody detects endogenous levels of Ribosomal Protein L39L protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	60S ribosomal protein L39-2;60S ribosomal protein L39-like;Large ribosomal subunit protein eL3;Ribosomal Protein L39L;RPL39L;RPL39L1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytosolic large ribosomal subunit,
Tissue Specificity	Testis specific.
Function	similarity:Belongs to the ribosomal protein L39e family.,tissue specificity:Testis specific.,
Background	This gene encodes a protein sharing high sequence similarity with ribosomal protein L39. Although the name of this gene has been referred to as 'ribosomal protein L39' in the public databases, its official name is 'ribosomal protein L39-like'. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome. [provided by RefSeq, Jul 2008],



**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemistry analysis of paraffin-embedded human brain, using RPL39L Antibody. The picture on the right is blocked with the synthesized peptide.