




SAA4 Polyclonal Antibody

Catalog No	YP-Ab-04184
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	SAA4
Immunogen	The antiserum was produced against synthesized peptide derived from human SAA4. AA range:81-130
Specificity	SAA4 Polyclonal Antibody detects endogenous levels of SAA4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C SAA;CSAA;SAA4;Serum amyloid A 4 protein;serum amyloid A4;Serum amyloid A-4 protein;Constitutively expressed serum amyloid A protein;C-SAA
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Expressed by the liver; secreted in plasma.
Function	function:Major acute phase reactant. Apolipoprotein of the HDL complex., induction:Constitutively expressed.,similarity:Belongs to the SAA family.,tissue specificity:Expressed by the liver; secreted in plasma.,
Background	function:Major acute phase reactant. Apolipoprotein of the HDL complex., induction:Constitutively expressed.,similarity:Belongs to the SAA family.,tissue specificity:Expressed by the liver; secreted in plasma.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma, using SAA4 Antibody. The picture on the right is blocked with the synthesized peptide.