



TAFI Polyclonal Antibody

Catalog No	YP-Ab-04236
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	CPB2
Immunogen	The antiserum was produced against synthesized peptide derived from human CPB2. AA range:321-370
Specificity	TAFI Polyclonal Antibody detects endogenous levels of TAFI protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SET;2PP2A;HLA-DR-associated protein II;I-2PP2A;Inhibitor of granzyme A-activated DNase;TAFI;CPB2;Carboxypeptidase B2;Carboxypeptidase U;CPU; Plasma carboxypeptidase B;pCPB;Thrombin-activable fibrinolysis inhibitor
Observed Band	48kD
Calculated Molecular Weight	48kD
Cell Pathway	Secreted.
Tissue Specificity	Plasma; synthesized in the liver.
Function	catalytic activity:Release of C-terminal Arg and Lys from a polypeptide., cofactor:Binds 1 zinc ion per subunit.,function:Cleaves C-terminal arginine or lysine residues from biologically active peptides such as kinins or anaphylatoxins in the circulation thereby regulating their activities., similarity:Belongs to the peptidase M14 family.,tissue specificity:Plasma; synthesized in the liver.,



Background

Carboxypeptidases are enzymes that hydrolyze C-terminal peptide bonds. The carboxypeptidase family includes metallo-, serine, and cysteine carboxypeptidases. According to their substrate specificity, these enzymes are referred to as carboxypeptidase A (cleaving aliphatic residues) or carboxypeptidase B (cleaving basic amino residues). The protein encoded by this gene is activated by trypsin and acts on carboxypeptidase B substrates. After thrombin activation, the mature protein downregulates fibrinolysis. Polymorphisms have been described for this gene and its promoter region. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jun 2013],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from HepG2 cells, using CPB2 Antibody. The lane on the right is blocked with the synthesized peptide.

Western blot analysis of the lysates from K562 cells using CPB2 antibody.