




Tctex5 Polyclonal Antibody

Catalog No	YP-Ab-04240
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	IHC;IF;ELISA
Gene Name	PPP1R11
Immunogen	The antiserum was produced against synthesized peptide derived from human PPP1R11. AA range:71-120
Specificity	Tctex5 Polyclonal Antibody detects endogenous levels of Tctex5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/40000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HCG V;HCGV;PPP1R11;PPP1R11-Specific;TCTE5;TCTEX5;Protein phosphatase 1 regulatory subunit 11;Hemochromatosis candidate gene V protein;Protein phosphatase inhibitor 3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	protein phosphatase type 1 complex,nucleus,cytoplasm,
Tissue Specificity	Widely expressed.
Function	function:Inhibitor of protein phosphatase 1.,tissue specificity:Widely expressed.,
Background	This gene encodes a specific inhibitor of protein phosphatase-1 (PP1) with a differential sensitivity toward the metal-independent and metal-dependent forms of PP1. The gene is located within the major histocompatibility complex class I region on chromosome 6. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using PPP1R11 Antibody. The picture on the right is blocked with the synthesized peptide.